blue dot[®]

B-BM03 Pro light

Automatic Digital

EN BLOOD PRESSURE MONITOR

User manual.....2

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CE 1639

EN : AUTOMATIC DIGITAL BLOOD PRESSURE MONITOR

Thank you for purchasing Automatic Digital Blood Pressure Monitor. It is mainly designed for general home use. Before using the device, please read this manual carefully to ensure proper and safe operation. Please take good care of the manual for future reference.

We glad to give you advice and support.

Table of Contents

1. Getting to know your device.....	03	5. Measuring blood pressure.....	10
2. Important notes.....	03	Attaching the cuff.....	10
Signs and symbols.....	03	Adopting the correct posture.....	10
Notes on use.....	04	Performing the blood pressure	
Instructions for storage and maintenance....	06	measurement.....	11
Instructions for repairs and disposal.....	06	Reading Memory.....	12
3. Device description.....	07	Deleting all memory values.....	12
4. Preparing for the measurement.....	08	Motion detector.....	12
Inserting the batteries.....	08	Evaluating results.....	12
Operation with the mains adapter.....	09	6. Abnormal phenomenon and Solutions.....	14
Display Setup and Memory Deletion.....	09	7. Maintenance and storage.....	14
Cuff Connection.....	09	8. Technical Specifications.....	15

Included in delivery

- Blood pressure monitor 1 pc
- Upper arm cuff size 22-42 cm 1 pc
- 1.5 V AA batteries 4 pcs
- Adapter 1 pc
- Storage bag 1 pc
- User manual 1 pc

1. Getting to know your device

Check that the device packaging has not been tampered and make sure that all contents are present. Before use, ensure that there is no visible damage to the device or accessories and that all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Services address. This new digital blood pressure monitor is an automatic measuring device to be used on the upper arm. It enable simple, accurate and fast measurement of blood pressure and pulse rate by use of the oscillometric method and people-orientated design.

The advantages of this product :

1. Talking Thai Language
2. LCD screen with backlight and 3 color display

3. Intel comfort
4. Irregular heartbeat detection
5. Motion detection
6. 60 Memories / person (2 person)
7. WHO Indicator
8. Adapter and Battery



2. Important notes

Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and the accessories:

	Attention
	Note Note on important information
	Observe the instructions for use
	Application part, type BF

	Direct current
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE
	Dispose of packaging in an environmentally friendly manner
Storage/Transport 	Permissible storage and transport temperature and humidity
Operating 	Permissible operating temperature and humidity
	Protect from moisture
	The CE labeling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products.

Notes on use

- In order to ensure comparable values, always measure your blood pressure at the same time of day.
- Before every measurement, relax for about five minutes.
- If you want to perform several measurements on the same person, wait five minutes between each measurement.
- Do not take a measurement within 30 minutes of eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.
- The measured values taken by you are for your information only – they are no substitute for a medical examination. Discuss the measured values with your doctor and never base any medical decisions on them (e.g. medicines and their administration).
- Using the blood pressure monitor outside your home environment or whilst on the move (e.g. whilst travelling in a car, ambulance or helicopter, or whilst undertaking physical activities such as playing sport) can influence the measurement accuracy and cause incorrect measurements.

- Do not use the blood pressure monitor on newborns or patients with pre-eclampsia. We recommend consulting a doctor before using the blood pressure monitor during pregnancy.
- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.
- This device is not intended for use by people (including children) with restricted physical, sensory or mental skills or a lack of experience and/or a lack of knowledge, unless they are supervised by a person who is responsible for their safety or are instructed by such a person in how to use the device. Supervise children around the device to ensure they do not play with it.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the device on people who have the specified upper arm measurement for the device.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, the blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions remove the cuff from the arm.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.
- Make sure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or intravascular or therapy, or an arteriovenous (AV) shunt.
- Do not use the cuff on people who have undergone a mastectomy.
- Do not place the cuff over wounds as this may cause further injury.
- Place the cuff on your upper arm only. Do not place the cuff on other parts of the body.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.

Instructions for storage and maintenance

- The blood pressure monitor is made from precision and electronic components. The accuracy of the measured values and service life of the device depend on its careful handling:
 - Protect the device from impacts, humidity, dirt, marked temperature fluctuations and direct sunlight.
 - Do not drop the device.
 - Do not use the device in the vicinity of strong electromagnetic fields and keep it away from radio systems or mobile telephones.
 - Only use the cuff included with the delivery or original replacement parts. Otherwise incorrect measured values will be recorded.
- We recommend that the batteries be removed if the device will not be used for a prolonged period of time.

Notes on handling batteries

- If your skin or eyes come into contact with battery fluid, rinse the affected areas with water and seek medical assistance.
-  **Choking hazard!** Small children may swallow and choke on batteries. Store the batteries out of the reach of small children.

- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect batteries from excessive heat.

 **Risk of explosion!** Never throw batteries into a fire.

- Do not charge or short-circuit batteries.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, split or crush the batteries.

Instructions for repairs and disposal

- Batteries do not belong in household waste. Please dispose of empty batteries at the collection points intended for this purpose.
- Do not open the device.
Failure to comply will result in voiding of the warranty.
- Do not repair or adjust the device yourself. Proper operation can no longer be guaranteed in this case.
- Repairs must only be carried out by Customer Services or authorised suppliers. Before making a claim, please check the batteries first and replace them if necessary.


- Please dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

Notes on electromagnetic compatibility


- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of accessories other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Failure to comply with the above can impair the performance of the device.

3. Device description

1. Cuff
2. Cuff line
3. Cuff connector
4. Display
5. Connection for cuff connector
6. Risk indicator (WHO standard)
7. Thai Voice Symbol
8. Intel comfort Symbol
9. Memory button 
10. START/STOP button 
11. Setting button 
12. AC adapter Support

Information on the display:


1. Systolic pressure (SYS) (mmHg)
2. Diastolic pressure (DIA) (mmHg)
3. Calculated pulse value
4. Time and date
5. Cardiac arrhythmia symbol (♥)
6. Average of the last 3 measurements
7. User memory (👤 👤)
8. Motion detector (👤)
9. Battery replacement symbol (🔋)
10. Off/On Voice (🔊)
11. Pulse symbol (♥)
12. Risk indicator (WHO standard)

4. Preparing the measurement

Inserting the batteries

Press the indicator ▲ on the battery cover and slide the cover off in the direction of the arrow. Insert 4 "AA" size batteries, align the polarities of the battery.

Note: If the low battery indicator appears on the display, the batteries must be replaced. This monitor can also be operated with an AC adapter (DC 6V).

WARNING Please purchase AC adapter with CE certification. Do not plug or unplug the power cord into the electrical outlet with wet hands.

If the low battery indicator is permanently displayed, you can no longer perform any measurements and must replace the batteries.

Once the batteries have been removed from the device, the time must be set again.

Battery disposal

- The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

Operation with the mains adapter

You can also operate this device with a mains adapter. When doing so, there must not be any batteries in the battery compartment.

- To prevent possible damage to the device, the blood pressure monitor must only be used with the mains adapter described here.
- Insert the mains adapter into the connection provided for this purpose on the blood pressure monitor. The mains adapter must only be connected to the mains voltage that is specified on the type plate.
- Then insert the mains plug of the mains adapter into the mains socket.
- After using the blood pressure monitor, unplug the mains adapter from the mains socket first and then disconnect it from the blood pressure monitor.

As soon as you unplug the mains adapter, the blood pressure monitor loses the date and time setting but the saved measurements are retained.

Display Setup and Memory Deletion

It is essential to set the date and time. Otherwise, you will not be able to save your measured values correctly with a date and time and access them again later.

- Before initial use and after each time you replace the battery: When inserting batteries into the device, you will be taken to the relevant menu automatically.
- If the batteries have already been inserted: Press and hold the settings button SET  on the device when switched off for approx. 3 seconds. Enter into the setup interface, the content of setup will be displayed by flashing. Press the SET button, it will enter into the switch of setup content:

"memory deletion -> talking function -> year -> month -> day -> hour -> minute."

- Press memory button  to revise setup value.
- Press ON/OFF button  to store the setup content.

Note: This monitor displays one unit mmHg.

Cuff Connection


Insert the air tube plug into the socket by revolving it.

5. Measuring blood pressure


Ensure the device is at room temperature before measuring. The measurement can be performed on the left or right arm.

Attaching the cuff


Place the cuff onto the bare upper arm. The circulation of the arm must not be hindered by tight clothing or similar.

The cuff must be placed on the upper arm so that the bottom edge is positioned 1 – 2 cm above the elbow and over the artery. The line should point to the centre of the palm.

Now tighten the free end of the cuff, but make sure that it is not too tight around the arm and close the hook-and-loop fastener. The cuff should be fastened so that two fingers can fit under the cuff.

Blood pressure may vary between the right and left arm, which may mean that the measured blood pressure values are different. Always perform the measurement on the same arm. If the values between the two arms are significantly different, please consult your doctor to determine which arm should be used for the measurement.


i Important: The unit may only be operated with the original cuff. The cuff is suitable for an arm circumference of 22 to 42 cm.


Adopting the correct posture

Sit upright with your back straight.

Take off thick clothes, wear thin clothes or bare your upper arm for measurement.

The arm cuff should be at the same level as your heart.


- Before the initial blood pressure measurement, make sure always to rest for about 5 minutes. Otherwise deviations can occur.
- Furthermore, if you want to take several measurements in succession, make sure always to wait for at least 3-5 minutes between the individual measurements.
- You can take the measurement while sitting or lying. Always make sure that the cuff is at heart level.
- To carry out a blood pressure measurement, make sure you are sitting comfortably with your arms and back leaning on something. Do not cross your legs. Place your feet flat on the ground.
- To avoid falsifying the measurement, it is important to remain still during the measurement and not to speak.

Performing the blood pressure measurement

1. Press SET button  to choose user A or user B.
2. Attaching the cuff
3. Press ON/OFF button  to power on, The pump starts to inflate the cuff.

4. After having reached the anticipated value of inflation, the monitor starts to deflate automatically with a constant speed to detect your blood pressure and pulse.

5. When the measurement is completed, your blood pressure and pulse rate will be displayed.

Example:

SYS : 128 mmHg
DIA : 78 mmHg
Pulse : 72 bpm

The result will be stored automatically.

The storage series number is "0", doesn't have memory.

The previous measurement serial number is "1,2,3,...,60"

6. To shut off the monitor, press ON/OFF button. 

 **Note:** Once you feel uncomfortable with inflation or the inflation becomes ceaseless, please loosen the cuff and pull out the air tube immediately. Don't talk, move your arm or body during the measurement.

The monitor will be automatically shut off in 2 minutes without operation. When the monitor can't measure the value correctly, the screen will display "E" icon. Press ON/OFF button to shut off the monitor and wait for 4 to 5 minutes to start a new measurement. At the same time, follow the following suggestions during the measurement to avoid "E" icon display again.

1. Keep quiet and correct posture.
2. Fit the cuff correctly.
3. Connect the air tube correctly.
4. Avoid the electromagnetic interference.

Reading Memory

Press memory button, the average value of latest three times measurement will be displayed with  symbol.

Press memory button  to read memory.

Press memory button or SET button  to cycle switch.

Note: The monitor can memorize 60 data.

The monitor will delete the earliest data automatically when the number of data is beyond 60.

Deleting all memory values

Press SET button  for about 3 seconds to enter the setting mode. All memory values show flashing status. Press Memory  to select Delete. "0" will delete all memory values and press the On/Off button  to confirm.

Motion detector

The icon will show on display when device detect movement during measurement.

Evaluating results

Cardiac arrhythmia:

This device can identify potential disruptions of the heart rhythm when measuring and if necessary, indicates this after the measurement with the symbol . This can be an indicator for arrhythmia.

Arrhythmia is an illness in which the heart rhythm is abnormal because of flaws in the bioelectrical system that regulates the heartbeat. The symptoms (skipped or premature heart beats, pulse being slow or too fast) can be caused by factors such as heart disease, age, physical make-up, excess stimulants, stress or lack of sleep. Arrhythmia can only be determined through an examination by your doctor. If the symbol  is shown on the display after the measurement has been taken, repeat the measurement.

Please ensure that you rest for 3-5 minutes beforehand and do not speak or move during the measurement. If the symbol  appears frequently, please consult your doctor. Self diagnosis and treatment based on the measurements can be dangerous. Always follow your GP's instructions.

Risk indicator:

The measurements can be classified and evaluated in accordance with the following table. However, these standard values serve only as a general guideline, as the individual blood pressure varies in different people and different age groups etc. It is important to consult your doctor regularly for advice. Your doctor will tell you your individual values for normal blood pressure as well as the value above which your blood pressure is classified as dangerous.

Table for Classifying Blood Pressure Values (mmHg): WHO

Blood Pressure Classification	Risk indicator	SYS mmHg	DIA mmHg	Action
Severe Hypertension	Red	≥180	≥110	seek medical attention
Hypertension	Yellow	160-179	100-109	seek medical attention
Mild Hypertension	Yellow	140-159	90-99	regular monitoring by doctor
High-normal	Green	130-139	85-89	regular monitoring by doctor
Normal		120-129	80-84	self-monitoring
Optimal	Green	<120	< 80	self-monitoring

Life Guide for Hypertension Patients

• High blood pressure (Hypertension);
A blood pressure of 140/90 or higher is considered high blood pressure. High blood pressure (HBP) is a serious condition that can lead to coronary heart disease, heart failure, stroke, kidney failure, and other health problems. Blood pressure tends to rise with age. Following a healthy lifestyle helps some people delay or prevent this rise in blood pressure. The following measures can be regarded as a life guide for hypertension patients:

- 1) Consult your doctor about your blood pressure in daily life.
- 2) Maintain a normal weight.
- 3) Avoid excessive consumption of common salt. Avoid high cholesterol foods and fatty foods. Eat more fruits, vegetables and high-fiber foods.
- 4) Avoid smoking, alcohol and reviver.
- 5) Moderated exercise.
- 6) Relax and reduce stress. Stress can increase the blood pressure suddenly. Relax and reduce stress is a good habit for your health.

6. Abnormal phenomenon and Solutions

Error Indicators

Symbol	Cause	Solutions
E1	Attaching the cuff in wrong position	Please check the correct attaching the cuff see page 10.
	Air tube plug is not correctly connected	Check the connection of the air tube and connect it properly if necessary. Restart measuring.
EE	Systolic value above 210 mmHg	Please don't move and speak while measuring or contact your distributor
	The inflation is under 30 mmHg	
E4	The monitor or arm is moved during the measurement	Keep your body still, never move the monitor. Restart measuring
E5	The inflation is above 300 mmHg	Please don't move and speak while measuring or contact your distributor

Troubleshooting

Abnormal Phenomenon	Potential cause	Solutions
Power on doesn't work	Low battery or wrong polarities of batteries aligned	Replace old batteries with new ones; Align the polarities of the battery correctly.
No inflation	Air tube plug is not correctly connected	Check the connection of the air tube and connect properly if necessary.
Power off when inflation	Low battery	Replace old batteries with new ones
Abnormal value	Cuff is not fitted on the arm correctly	Fit the cuff correctly see page 10.
	Talk during the measurement	Keep quiet during measurement

Abnormal value	Roll-up sleeve pressed on the arm	Take off the clothes on the arm, restart measuring
	Nervous or excited	Keep quiet, take a deep breath to calm down
	Wrong measurement	Adjust the posture see page 11

7. Maintenance and storage

- Do not subject the monitor and package to shocks, such as dropping it onto the floor.
- Do not submerge the device or any of the components in water. Do not subject the monitor to extreme temperatures, dust, humidity, corrosive gas and direct sunlight, do not use it in an oxygen rich environment.
- Use the cuff carefully. The cuff contains a sensitive gasbag and should not be forcibly bended.
- Clean the monitor with soft and dry cloth. Don't use gasoline, thinners or similar solvents.
- Carefully remove spots on the cuff with a damp cloth and soapsuds. It should not be washed with detergent powder and other detergents.
- Use only approved parts and accessories. Not approved parts and accessories may damage the unit.
- Store the device and the components in a clean, safe location.
- About batteries: If the monitor is not to be used for a long time, remove the batteries to prevent battery leakage.

8. Technical Specifications

Model No.	B-BM03 Pro light
System	Fully Automatic inflation and deflation
Measuring Range	Cuff Pressure: 0-290 mmHg Pulse: 40~180 bpm
Accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$ of the reading
Display	Illuminated LCD digital display 65 x 90 mm
Memory	2 x 60 memory spaces
Dimension	110 x 150 x 70 mm
Weight	310 g (without batteries)
Cuff size	22-42 cm
Operating Condition	Temperature: 5°C~40°C Humidity: 15% RH~85% RH Air Pressure: 860 hPa~1060 hPa
Storage Condition	Temperature: -20°C~70°C Humidity: 10% RH~95% RH Air Pressure: 860 hPa~1060 hPa
Power source	4X1.5V "AA" batteries or AC adapter DC6V

Guidance and Manufacturer's Declaration

Important information regarding Electro Magnetic Compatibility (EMC)

This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.

- * Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
- * Caution: this machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Guidance and manufacturer's declaration - electromagnetic emission		
The device is intended for use in the electromagnetic environment specified below. The customer of the user of the device should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable (Battery operated device)	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable (Battery operated device)	

Guidance and manufacturer's declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor is covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not applicable (Battery operated device)	Mains power quality should be that of a typical commercial or hospital environment.

Surge IEC 61000-4-5	± 1 kV /line(s) to line(s) ± 2 kV /line(s) to earth	Not applicable (Battery operated device)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T for 5 sec (>95% dip in U_T)	Not applicable (Battery operated device)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital

Guidance and manufacturer's declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF JEC 61000-4-6</p>	<p>3 V_{rms} 150 kHz to 80 MHz</p>	<p>Not applicable (Battery operated device) 3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> <p>$d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2.5 GHz</p>
<p>Radiated RF IEC 61000-4-3</p>	<p>3 V/m 80 MHz to 2.5 GHz</p>		<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: </p>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rate maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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
Guidance and manufacture's declaration - electromagnetic emission		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable (Battery operated device)	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable (Battery operated device)	

Guidance and manufacture's declaration - electromagnetic Immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor is covered with synthetic material, the relative humidity should be at least 30%.

Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable (Battery operated device)	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable (Battery operated device)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$<5\% U_T$ ($>95\%$ dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $<5\% U_T$ for 5 sec ($>95\%$ dip in U_T)	Not applicable (Battery operated device)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital

Guidance and manufacture's declaration - electromagnetic Immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	$3 V_{rms}$ 150 kHz to 80 MHz	Not applicable (Battery operated device)	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ² should be less than the compliance level in each frequency range. ³ Interference may occur in the vicinity of equipment marked with the following symbol: 

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LIFEBOX ASIA.Co.,Ltd