blue●dot®

PERSONAL SCALE

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USER MANUAL: EN

Dear user,

Your new electronic personal scale has been designed to accurately indicate your weight gain or loss. Read these instructions for use carefully and keep them for later use.

SPECIFICATIONS

- Product size 350 X 300 X 24 mm
- Equipped with a High Precision "Strain-Gauge" Sensor
- 6 mm tempered safety glass platform
- Capacity 180 kg / 400 lb
- Graduation 0.1 kg / 0.2 lb
- LCD Display size 96 X 82 mm with backlight
- Auto-on/Auto-off function
- Tare function
- Overload Indication
- Low battery Indication
- Power 3XAAA batteries

OPERATION

- STEP 1. Put the scale on a hard & flat surface (avoid carpet or soft surface).
- STEP 2. Gently step on the scale. Put your feet evenly on the scale without moving and wait while the scale computes your weight. Your correct weight will then be displayed as shown in the picture:

- **Tare function:**

5 seconds after the first measurement, the screen will reset to 0.0 kg and display "Tare". Within the next 30 seconds, step on the scale again while carrying an additional load, such as a pet. The screen will then show the weight of the added load.

AUTOMATIC SWITCH OFF

When you step off the scale, it will automatically switch off. (if no weighing is carried out while the zero is being displayed, the scale switches off automatically)

OVERLOAD INDICATOR

When the scale is over loaded (over 180 kg / 400 lb), the overload indicator will appear on the screen.

LOW BATTERY INDICATOR

When the battery is running out, the low battery indicator will appear on the screen. Please open the battery cover and replace new batteries.

ADVICE FOR USE AND CARE

- Place the scale on a flat surface. It will ensure the greatest accuracy.
- Always weigh yourself on the same scale placed on the same floor surface. Uneven floors may affect readings.
- Do not allow children to use the scale without adult supervision.
- You can clean the scale with a damp cloth and a little washing up liquid if required, but do not immerse the scale in water or rinse it under running water as this can damage the electronics.

- Always insert the batteries correctly and remove the batteries if the scale is not being used for a long time.
- Treat your scale with care. Do not drop it or jump on it.

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Product Owner: Lifebox Asia Company Limited